

### How to Frisk

- Always frisk slowly, at a rate of ~2 cm/s and
- Never touch the detector directly, always leave ~1 cm between the object you are frisking and the detector head.
- When frisking out, frisk both bottom of feet and hands. Only pick up frisker after frisking that hand.
- Be sure to **check any personal items** you are carrying (books, clipboards, backpacks, etc...) for contamination.

### Definitions

Absorbed Dose: The energy imparted per ions or charged particles.

unit mass or irradiated material. Measured in **rad**, where 1 rad equals 0.01 Joules/kg of absorbing material. The SI unit is the Gray (Gy).

Note 1 Gy = 100 rad =1J/kg

**Exposure:** A measure of the ionization produced in *AIR* by X or Gamma radiation. Measured in **Roentgen (R)**, where  $1R = 2.58 \times 10^4$  Coulombs per kilogram of dry air at STP.

<u>Gamma Radiation (y</u>): Short wavelength electromagnetic radiation of nuclear origin.

#### VERY high penetration ability.

<u>Ionizing Radiation</u>: Electromagnetic (X ray and gamma) or particulate (alpha, beta) radiation capable of producing **<u>Radiation</u>**: Energy transmitted as electromagnetic waves or particles from a source.

**Radioactivity:** The property of certain nuclides of spontaneously emitting particles or gamma radiations or emitting X-rays following orbital electron capture. Measured units are Curies (Ci) or SI unit, Becquerel's (Bq). 1 Bq = 1 disintegration per second. 1 Ci =  $3.7x10^{10}$  Bq.

RAM: Radioactive Material

# Rule for Irradiator Use

After completing the security entrance requirements:

- As the user or the escort, you must fill in the user log sheet. Be sure to fill in all fields and do NOT use "ditto" or check marks.
- Check the irradiation timer to ensure the irradiator is not currently in use.
- Survey the Chamber door prior to opening it.
- Do not remove someone else's sample unless authorized to do so.
- Answer the room phone when/if it rings.
- Wear your dosimeter.
- Clean compartment after use, remove all generated trash from the room.

### **Decay Correction Equation**

D(t) = dose rate at desired time  $D_o$  = reference dose rate  $\lambda$  = Ln(2)/ $T_{1/2}$   $T_{1/2}$  = half life of isotope t = time lapse in years

\* \*This formula is used to generate the dose rate table posted for each irradiator.


$$D(t) = D_o e^{-\lambda t}$$

#### Performing a Survey

Before and after working with the irradiator, always perform a survey of the chamber door before opening.

- Check calibration date on detector
- Check battery
- Check response using check source mounted on instrument.
- Measure the background exposure rate in an area known to have low radiation.
- Scan around the chamber door by holding detector face ~1 cm away from the surface moving slowly (about 2 inches/s).

### Notify RPP if unusually high exposure levels are detected



## 9 Traits of a Positive Safety Culture

- 1. Leadership Safety Values and Actions
- 2. Problem Identification
- 3. Personal Accountability
- 4. Work Processes
- 5. Continuous Learning

- 6. Environment for Raising Concerns
- 7. Effective Safety Communication
- 8. Respectful Work Environment
- 9. Questioning Attitude